

Product Name: Recombinant Human FGGY Protein
Catalog #: PEH9635



Summary

Name	FGGY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96C11
Host	E.coli
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YYVGVDVGTG SVRAALVDQS GVLLAFADQP IKNWEPQFNH HEQSSEDIWA ACCVVTKKV V QGIDLNQIRG LGFDATCSLV VLDKQFHPLP VNQEGDSHRN VIMWLDHRAV SQVNRINETK HSVLQYVGGV MSVEMQAPKL LWLKENLREI CWDKAGHFFD LPDFLSWKAT GVTARSLCSL VCKWTYSAEK GWDDSFWKMI GLEDFVADNY SKIGNQVLPP GASLGNGLTP EAARDLGLLP GIAVAASLID AHAGGLGVIG ADVRGHGLIC EGQPVTSRLA VICGTSSSCHM GISKDPIFVP GVWGPYFSAM VPGFWLNEGG QSVTGKLIDH MVQGHAAPFE LQVKATARCQ SIYAYLNSHL DLIKKAQPVG FLTVDLHVWP DFHGNRSPLA DLTGKGMVTG LKLSQDLDDL AILYLATVQA IALGTRFIE AMEAAGHSIS TLFLCGGLSK NPLFVQMHAD ITGMPVLSQ EVESVLVGAA VLGACASGDF ASVQEAMAKM SKVGKVVFPRLQDKKYYDKK YQVFLKLVEH QKEYLAIMND

image



Background

Product Name: Recombinant Human FGGY Protein
Catalog #: PEH9635



Alternative Names

FGGY carbohydrate kinase domain-containing protein (EC 2.7.1.-)

Note

For Research Use Only , Not for Diagnostic Use