


Product Name: Recombinant Human FHAD1 Protein
Catalog #: PEH9636



Summary

Name	FHAD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	B1AJZ9
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKETLEEERK RMQELESLLA QQKKALAKSI TQEKNRVKEA LEEEQTRVQE LEERLARQKE VLESSIAHEK RKAKEALESE KRKVQDLENH LTQQKEISES NIAYEKRKAK EAMEKEKKKV QDLENRLTKQ KEELELKEQK EDVLNKNLSD ALAMVEETQK TKATESLKAE SLALKLNETL AELETTKTKM IMVEERLILQ QKMKALQDE

image



Background

Alternative Names Forkhead-associated domain-containing protein 1 (FHA domain-containing protein 1)

Note

Product Name: Recombinant Human FHAD1 Protein
Catalog #: PEH9636



For Research Use Only , Not for Diagnostic Use