


Product Name: Recombinant Human FAHD2B Protein
Catalog #: PEH9638



Summary

Name	FAH2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6P2I3
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QFRAPHLVGP HLGLETGNGG GVINLNAFDP TLPKTMTQFL EQGEATLSVA RRALAAQLPV LPWSEVTFLA PVTWPKVVC VGMNYVDHCK EQNVPVPKEP IIFSKFASSI VGPYDEVVLP PQSQEVDWEV ELAVVIGKKG KHIKATDAMA HVAGFTVAHD VSARDWLTRR NGKQWLLGKT FDTFCPLGPA LVTKDSVADP HNLKICCRVN GEVVQSSNTN QMVFKTEDLI AWVSQFVTFY PGDVILTGTP PGVGVFRKPP VFLKKGDEVQ CEIEELGVII

image



Background

Alternative Names Fumarylacetoacetate hydrolase domain-containing protein 2B (EC 3.-.-.-)

Product Name: Recombinant Human FAHD2B Protein
Catalog #: PEH9638



Note

For Research Use Only , Not for Diagnostic Use