


Product Name: Recombinant Human FYCO1 Protein
Catalog #: PEH9639



Summary

Name	FYCO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BQS8
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGQQGVGPPT DNEARELAAQ LALSQAQLEV HQGEVQRLQA QVVDLQAKMR AALDDQDKVQ SQLSMAEAVL REHKTIVQQL KEQNEALNRA HVQELLQCSE REGALQEERA DEAQREEEL RALQEELSQA KCSSEEAQLE HAELQEQLHR ANTDTAELGI QVCALTVEKE RVEEALACAV QELQDAKEAA SRERGLERQ VAGLQQEKES LQEKLKAAKA AAGSLPGLQA QLAQAEQRAQ SLQEAAHQEL NTLKFQLSAE IMDYQSRLKN AGEECKSLRG QLEEQGRQLQ AAEAVEKLN ATQADMGEKL SCTSNHLAEC QAAMLRKDKE GAALREDLER TQKELEKATT KIQEYYNKLC QEVTNRERND QKMLADLDDL NRTKKYLEER LIELLRDKDA LWQKSDALEF QQKLSAERW LGDTEANHCL

image



Background

Product Name: Recombinant Human FYCO1 Protein
Catalog #: PEH9639



Alternative Names

FYVE and coiled-coil domain-containing protein 1 (Zinc finger FYVE domain-containing protein 7)

Note

For Research Use Only , Not for Diagnostic Use