


Product Name: Recombinant Human GPRIN2 Protein
Catalog #: PEH9644



Summary

Name	GRIN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60269
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGREQRPELR KTASSTVWQA QLGEASTRPQ APEEEGNPPE SMKPARASGP KARPSAGGHW WSSTVGNVST MGGSDLCRLR APSAAAMQRS HSDLVRSTQM RGHSGARKAS LSCSALGSSP VHRAQLQPGG TSGQGGQAPA GLERDLAPED ETSNSAWMLG ASQLSVPLD LGDTTAHSS

image



Background

Alternative Names G protein-regulated inducer of neurite outgrowth 2 (GRIN2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GPRIN2 Protein
Catalog #: PEH9644

