


Product Name: Recombinant Human GSE1 Protein
Catalog #: PEH9649



Summary

Name	GSE1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14687
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SHSSLAALHS ERMSGLSAER LQMDEELRRE RERERERERE READREREKE REREREKERE QEKEREREKE RERELERQRE QRAREKELLA AKALEPSFLP VAELHGLRGH ATEERGKPSE QLTPTRAEKL KDAGLQAPKP VQHPLHPVPT PHHTVPSLIS NHGIFSLPSS SAATALLIQR TNEEEKWLAR QRRLRQEKEKED RQSQVSEFRQ QVLEQHLDMG RPPVPAEAEH

image



Background

Alternative Names Genetic suppressor element 1

Note

Product Name: Recombinant Human GSE1 Protein
Catalog #: PEH9649



For Research Use Only , Not for Diagnostic Use