


Summary

Name	PGTB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P53611
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQLHRMNREE ILAFIKSCQH ECGGISASIG HDPHLLYTLS AVQILTLYDS INVIDVNKVV EYVKGLQKED GSFAGDIWGE IDTRFSFCAV ATLALLGKLD AINVEKAIEF VLSCMNFDDG FGCRPGSESH AGQIYCCTGF LAITSQLHQV NSDLLGWWLC ERQLPSGGLN GRPEKLPDVC YSWWVVLASLK IIGRLHWIDR EKLRNFILAC QDEETGGFAD RPGDMVDPFH TLFGIAGLSL LGEEQIKPVN PVFCMPPEEVL QRVNVQPELV

image



Background

Alternative Names	Geranylgeranyl transferase type-2 subunit beta (EC 2.5.1.60) (Geranylgeranyl transferase type II subunit beta) (GGTase-II-beta) (Rab geranylgeranyltransferase subunit beta) (Rab GG transferase beta) (Rab GGTase)
--------------------------	---

Product Name: Recombinant Human RABGGTB Protein
Catalog #: PEH9650



beta) (Rab geranylgeranyltransferase subunit beta) (Type II protein geranylgeranyltransferase subunit beta)

Note

For Research Use Only , Not for Diagnostic Use