


Product Name: Recombinant Human GAN Protein
Catalog #: PEH9651



Summary

Name	GAN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2C0
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VREIVKECSN IPLSQPQQGE AMLANFKPRG YSECIVTVGG EERSVRKPTA AMRCMCPLYD PNRQLWIELA PLSMPRINHG VLSAEGFLFV FGGQDENKQT LSSGEKYDPD ANTWTALPPM NEARHNFGIV EIDGMLYILG GEDGEKELIS MECYDIYSKT WTKQPDLMV RKIGCYAAMK KKIYAMGGGS YGKLFESVEC YDPRTQQWTA ICPLKERRFG AVACGVAMEL YVFGGVRSRE DAQGSEMVTC KSEFYHDEFK RWIYLNDQNL

image



Background

Alternative Names Gigaxonin (Kelch-like protein 16)

Product Name: Recombinant Human GAN Protein
Catalog #: PEH9651



Note

For Research Use Only , Not for Diagnostic Use