


Product Name: Recombinant Human GMEB2 Protein
Catalog #: PEH9652



Summary

Name	GMEB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKD1
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEVVVVVTPD TAVDGSGVEG VKTVLVTTNL APHGGDLTED NMETENAAAA AAAAFTASSQ LKEAVLVKMA EEGINLEAEI VYPITCGDSR ANLIWRKFVC PGINVKCVQY DEHVISPKF VHLAGKSTLK DWKRAIRMNG IMLRKIMDSG ELDFYQHDKV CSNTRSTKI DLSGARVSLS SPTSAEYIPL TPAAADVNGS PATITETCE DPGDWTAAIG DDTFTFWRGL KDAGLLDEVI QEFHQELVET MRGLQQRVQD PPLQLRDAVL LNNIVQNFGM LDLVKKVLAS HKCQMDSRE

image



Background

Alternative Names	Glucocorticoid modulatory element-binding protein 2 (GMEB-2) (DNA-binding protein p79PIF) (Parvovirus initiation factor p79) (PIF p79)
--------------------------	--

Product Name: Recombinant Human GMEB2 Protein
Catalog #: PEH9652



Note

For Research Use Only , Not for Diagnostic Use