


Product Name: Recombinant Human RAB3IL1 Protein
Catalog #: PEH9655



Summary

Name	R3GEF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TBN0
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GLPPPLAAMP VPWKSTDPCQ GHRESPGALV ETSAGEEAQG QEGPAAAQLD VLRRLRSSME IREKGSEFLK EELHRAQKEL KLKDEECERL SKVREQLEQE LEELTASLFE EAHKMVREAN MKQAASEKQL KEARGKIDML QAEVTALKTL

image



Background

Alternative Names	Guanine nucleotide exchange factor for Rab-3A (Rab-3A-interacting-like protein 1) (Rab3A-interacting-like protein 1) (Rabin3-like 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use