


Product Name: Recombinant Human HPS4 Protein
Catalog #: PEH9658



Summary

Name	HPS4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQG7
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGCILYKGLI VSTQLPPSLT AKVLLHRTAP QEQRLPTGED APQEHGAALP PNVQIIPVFV TKEEAI SLHE FPVEQMTRSL ASPAGLQDGS AQHHPKGGST SALKENATGH VESMAWTPD PTSPDEACPD GRKENGCLSG HDLESIRPAG LHNSARGEVL GLSSSLGKEL VFLQEELDLS

image



Background

Alternative Names Hermansky-Pudlak syndrome 4 protein (Light-ear protein homolog)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HPS4 Protein
Catalog #: PEH9658

