


Product Name: Recombinant Human HPS6 Protein
Catalog #: PEH9659



Summary

Name	HPS6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86YV9
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALELELLSS GRPKAVLQAV GQLVQKEQWD RALDAGLALG PSSPLLRSEI FKLLLAFAQ HRRLDAHLPL LCRLCPPELA PAELLLLRT YLPDEVPPT PFPEPGAEP LTVGLLKALL EQTGAQGWLS GVLSPYEDI LWDPTPPPT

image



Background

Alternative Names Hermansky-Pudlak syndrome 6 protein (Ruby-eye protein homolog) (Ru)

Note

For Research Use Only , Not for Diagnostic Use