


Product Name: Recombinant Human HSPBP1 Protein
Catalog #: PEH9661



Summary

Name	HPBP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9NZL4
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSAGGSGNSR PPRNLQGLLQ MAITAGSEEP DPPPEPMSEE RRQWLQEAMS AAFRGQREEV EQMKSLRVL SQPMPPTAGE AEQAADQQR

image



Background

Alternative Names Hsp70-binding protein 1 (HspBP1) (Heat shock protein-binding protein 1) (Hsp70-binding protein 2) (HspBP2) (Hsp70-interacting protein 1) (Hsp70-interacting protein 2)

Note

For Research Use Only , Not for Diagnostic Use