


## Summary

<b>Name</b>	HSDL2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q6YN16
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	46 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	CTVFITGASR GIGKAIALKA AKDGANIVIA AKTAQPHPKL LGTIYTAEE IEAVGGKALP CIVDVRDEQQ ISAAVEKAIK KFGGIDILVN NASAISLTNT LDTPTKRLDL MMNVNTRGTY LASKACIPYL KSKVAHILN ISPPLNLPV WFKQHCAITI AKYGMSMYVL GMAEEFKGEI AVNALWPKTA IHTAAMDMLG GPGIESQCRK VDIIADAAYS IFQKPKSFTG NFVIDENILK EEGIENFDVY AIKPGHPLQP DFFLDEYPEA VSKKVESTGA VPEFKEEKLQ LQPKPRSGAV EETFRIVKDS LSDDVVKATQ AIYLFELSGE DGGTWFLDLK SKGGNVGYGE PSDQADVVMMS MTTDDFVKMF SGCLKPTMAF MSGCLKIKGN MALAIKLEKL

## image



## Background

**Alternative Names** Hydroxysteroid dehydrogenase-like protein 2 (EC 1.-.-.-) (Short chain

**Product Name: Recombinant Human HSDL2 Protein**  
**Catalog #: PEH9662**



---

dehydrogenase/reductase family 13C member 1)

**Note**

For Research Use Only , Not for Diagnostic Use