


Product Name: Recombinant Human MINDY4B Protein
Catalog #: PEH9663



Summary

Name	MIIY4B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A8MYZ0
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSEQLDLEEI SRKISFLDKW REIFSYHRLG TNNSTPQNHE GNHTSADENE DGTGLSQPKG QGHLPSSGLC SIPNPSIISS

image



Background

Alternative Names Inactive ubiquitin carboxyl-terminal hydrolase MINDY-4B (Protein FAM188B2)

Note

For Research Use Only , Not for Diagnostic Use