


Product Name: Recombinant Human IPPK Protein
Catalog #: PEH9664



Summary

Name	IPPK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H8X2
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LYNRVERYLE EFPEERKTLQ IDGPYDEAFY QKLLDLSTED DGTVAFALTK VQQYRVAMTA KDCSIMIALS PCLQDASSDQ RPVVPSSRSR FAFSVSVLDL DLKPYESIPH QYKLDGKIVN YYSKTVRAKD NAVMSTRFKE

image



Background

Alternative Names	Inositol-pentakisphosphate 2-kinase (EC 2.7.1.158) (IPK1 homolog) (Inositol-1,3,4,5,6-pentakisphosphate 2-kinase) (Ins(1,3,4,5,6)P5 2-kinase) (InsP5 2-kinase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IPPK Protein
Catalog #: PEH9664

