


Product Name: Recombinant Human INTS3 Protein
Catalog #: PEH9665



Summary

Name	INT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q68E01
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPVSMEMDNH MSDKDESCYD NAEAAFSDDE EDLNSKGKKR EFRFHPIKET VVEEPVDITP YLDQLDESLR DKVLQLQKGS DTEAQCEVMQ EIVDQVLEED FDSEQLSVLA SCLQELFAH FRGEVLPEEI TEESLEESVG KPLYLIFRNL CQMQEDNSSF SLLDLLSEL YQKQPKIGYH LLYLRASKA AAGKMNLYES FAQATQLGDL HTCLMMDMKA

image



Background

Alternative Names	Integrator complex subunit 3 (Int3) (SOSS complex subunit A) (Sensor of single-strand DNA complex subunit A) (SOSS-A) (Sensor of ssDNA subunit A)
--------------------------	---

Note

Product Name: Recombinant Human INTS3 Protein
Catalog #: PEH9665



For Research Use Only , Not for Diagnostic Use