


Product Name: Recombinant Human IFFO1 Protein
Catalog #: PEH9668



Summary

Name	IFFO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q0D2I5
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LYNVLAKVKR ERDEYKRRWE EEYTVRIQLQ DRVNELQEEA QEADACQEEL ALKVEQLKAE LVVFKGLMSN NLSELDTKIQ EKAMKVDMDI CRRIDITAKL CDVAQQRNCE DMIQMFQVPS MGGRKRERKA AVEEDTSLSE SEGPRQPDGD EEESTALSIN EEMQRMLNQL REYDFEDDCD SLTWEETEET LLLWEDFSGY AMAAAEAQGE QEDSLEKVIK

image



Background

Alternative Names Intermediate filament family orphan 1 (Tumor antigen HOM-TES-103)

Note

Product Name: Recombinant Human IFFO1 Protein
Catalog #: PEH9668



For Research Use Only , Not for Diagnostic Use