


Product Name: Recombinant Human IFT122 Protein
Catalog #: PEH9669



Summary

Name	IF122_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HBG6
Host	E.coli
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VKKIAIYRNR LAIQLPEKIL IYELYSEDLS DMHYRVKEKI IKKFECNLLV VCANHIILCQ EKRLQCLSFS GVKEREWQME SLIRYIKVIG GPPGREGLLV GLKNGQILKI FVDNLFAIVL LKQATAVRCL DMSASRKKLA VVDENDTCLV YDIDTKELLF QEPNANSVAW NTQCEDMLCF SGGGYLNIKA STFPVHRQKL QGFVVGYNCS KIFCLHVFSI SAVEVPQSAP MYQYLDRKLF KEAYQIACLG VTDTDWRELA MEALEGLDFE TAKKAFIRVQ DLRYLELISS IEERKKRGET NNDLFLADV SYQGKFHEAA KLYKRSGHEN LALEMYTDLG MFEYAKDFLG SGDPKETKML ITKQADWARN IKEPKAAVEM YISAGEHVKA IEICGDHGWV DMLIDIARKL DKAEREPLLL CATYLLKLLS PGYAAETYLK MGDLSLVQL HVETQRWDEA FALGEKHPEF KDDIYMPYAQ WLAENDRFEE AQKAFHKAGR QREAVQVLEQ

image



Background

Product Name: Recombinant Human IFT122 Protein
Catalog #: PEH9669



Alternative Names

Intraflagellar transport protein 122 homolog (WD repeat-containing protein 10) (WD repeat-containing protein 140)

Note

For Research Use Only , Not for Diagnostic Use