


Product Name: Recombinant Human IFT122 Protein
Catalog #: PEH9670



Summary

Name	IF122_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HBG6
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LHLVEFYLEE GITDEEAI SL IDLEVLRPKR DDRQLEIANN SSQILRLVET KDSIGDEDPF TAKLSFEQGG SEFVPVVVSR LVLRSMSSRRD

image



Background

Alternative Names Intraflagellar transport protein 122 homolog (WD repeat-containing protein 10) (WD repeat-containing protein 140)

Note

For Research Use Only , Not for Diagnostic Use