


**Product Name: Recombinant Human TTC26 Protein**  
**Catalog #: PEH9672**



## Summary

<b>Name</b>	IFT56_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	A0AVF1
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	47 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	TEDQLSLASI HYMRSHYQEA IDIYKRILLD NREYLALNVY VALCYKLDY YDVSQEVAVLAV YLQQIPDSTI ALNLKACNHF RLYNGRAAEA ELKSLMDNAS SSFEFANELI RHNLVVFRGG EGALQVLPLP VDVIPPEARLN LVIYYLRQDD VQEAYNLIKD LEPTTPQEYI LKGVVNAALG QEMGSRDHMK IAQQFFQLVG GSASECDTIP GRQCMASCFE LLKQFDDVLI YLNSFKSYFY NDDIFNFNVA QAKAATGNTS EGEEAFLLIQ SEKMKNYIY LSWLARCIM NKKPRLAWEL YMKMETSSES FSLQLIAND CYKMGQFYYS AKAFDVLRL DPNPEYWEGK RGACVGIFQM IAGREPKEK LREVLHLLRS TGNTQVEYMI

## image



## Background

**Alternative Names**      Intraflagellar transport protein 56 (Tetratricopeptide repeat protein 26) (TPR)

**Product Name: Recombinant Human TTC26 Protein**  
**Catalog #: PEH9672**



---

repeat protein 26)

**Note**

For Research Use Only , Not for Diagnostic Use