


Product Name: Recombinant Human JAKMIP2 Protein
Catalog #: PEH9675



Summary

Name	JKIP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96AA8
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AQEQLQAEVL RYKAKIEDLE ATLAQKGQDS HWVEDKQLFI KRNQELLEKI EKQEAENHRL QQELQDARDQ NELLEFRNLE LEERERRSPP FNLQIHPFSD GVSALQIYCM KEGVKDVNIP DLIKQLDILG DNGNLRNEEQ VAIQASTVL SLAEKWIQQI EGAEAAALHQQ MMELESDMEQ FCKIKGYLEE ELDYRKQALD QAYMRIQELE ATLYNALQQE TVIKFGELLS EKQQEELRTA VEKLRRQMLR KSREYDCQIL QERMELLQQA

image



Background

Alternative Names	Janus kinase and microtubule-interacting protein 2 (CTCL tumor antigen HD-CL-04) (Neuroendocrine long coiled-coil protein 1)
--------------------------	--

Product Name: Recombinant Human JAKMIP2 Protein
Catalog #: PEH9675



Note

For Research Use Only , Not for Diagnostic Use