


Product Name: Recombinant Human KLHL12 Protein
Catalog #: PEH9679



Summary

Name	KLH12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q53G59
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTNTHAKSIL NSMNSLRKSN TLCDVTLRVE QKDFPAHRIV LAACSDYFCA MFTSESEKQ KPYVDIQGLT ASTMEILLDF VYTETVHVTV ENVQELLPAA CLLQLKGVKQ ACCEFLESQ L DPSNCLGIRD FAETHNCVDL MQAAEVFSQK HFPEVVQHEE FILLSQGEVE KLIKDEIQV DSEEPVFEAV INWVKHAKKE REESLPNLLQ YVRMPLTPR YITDVIDAEP FIRCSLQCRD

image



Background

Alternative Names Kelch-like protein 12 (CUL3-interacting protein 1) (DKIR homolog) (hDKIR)

Note

Product Name: Recombinant Human KLHL12 Protein
Catalog #: PEH9679



For Research Use Only , Not for Diagnostic Use