


Product Name: Recombinant Human KLHL31 Protein
Catalog #: PEH9681



Summary

Name	KLH31_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H511
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVEDSPLNKL NALNGLLEGG NGLSCISSEL TDASYGPNLL EGLSKMRQEN FLCDLVIGTK TKSFDVHKSV MASCSEYFYN ILKKDPSIQR VDLNDISPLG LATVIAYAYT GKLTLSLYTI GSIISAAVYL QIHTLVKMCS DFLIREMSVE NCMYVVNIAE TYSLKNAKAA AQKFIRDNFL EFAESDQFMK LTFEQINELL IDDDLQLPSE IVAFQIAMKW LEFDQKRVKY AADLLSNIRF GTISAQDLVN YVQSVPRMMQ DADCHRLDVD

image



Background

Alternative Names	Kelch-like protein 31 (BTB and kelch domain-containing protein 6) (Kelch repeat and BTB domain-containing protein 1) (Kelch-like protein KLHL)
--------------------------	--

Product Name: Recombinant Human KLHL31 Protein
Catalog #: PEH9681



Note

For Research Use Only , Not for Diagnostic Use