


Product Name: Recombinant Human KLHL8 Protein
Catalog #: PEH9682



Summary

Name	KLHL8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2G9
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVRAFAESHN RIDLMDMADQ YACDHFTEVV ECEDFVSVSP QLHLKLLSS DLNIENEKQV YNAAIKWLLA NPQHHSKWLD ETLAQVRLPL LPVDFLMGVV AKEQIVKQNL KCRDLLDEAR NYHLHLSSRA VPDFEYSIRT TPRKHTAGVL FCVGGRGGSG DPFRSIECYS INKNSWFFGP EMNSRRRHVG VISVEGKVYA VGGHDGNEHL GSMEMFDPLT NKWMMKASMN TKRRGIALAS LGGPIYAIGG LDDNTCFNDV ERYDIESDQW STVAPMNTPR GGVGSVALVN HUYAVGGNDG MASLSSVERY DPHLDKWIEV KEMGQRRAGN GVSKLHGCLY VVGGFDDNSP LSSVERYDPR SNKWDYVAAL

image



Background

Alternative Names Kelch-like protein 8

Product Name: Recombinant Human KLHL8 Protein
Catalog #: PEH9682



Note

For Research Use Only , Not for Diagnostic Use