


Product Name: Recombinant Human FN3KRP Protein
Catalog #: PEH9683



Summary

Name	KT3K_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HA64
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AILKTNTVKV PKPIKVLDP GGSVLVMEH MDMRHLSSHA AKLGAQLADL HLDNKKLGEM RLKEAGTVGR GGGQEERPFV ARFGFDVVC CGYLPQVNDW QEDWVVFYAR QRIQPQMDMV EKESGDREAL QLWSALQLKI PDLFRDLEII PALLHGDLWG GNVAEDSSGP

image



Background

Alternative Names	Ketosamine-3-kinase (EC 2.7.1.172) (Fructosamine-3-kinase-related protein) (FN3K-RP) (FN3K-related protein) (Protein-psicosamine 3-kinase FN3KRP) (EC 2.7.1.-)
--------------------------	--

Note

Product Name: Recombinant Human FN3KRP Protein
Catalog #: PEH9683



For Research Use Only , Not for Diagnostic Use