


Product Name: Recombinant Human KPTN Protein
Catalog #: PEH9684



Summary

Name	KPTN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y664
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PCPLREDSFT RFSSQSNVYG LAGGAGGRGE LLAATLKGKV LGFRYQDLRQ KIRPVAKELQ FNYIPVDAEI VSIDTFNKSP PKRGLVVGIT FIKDSGDKGS PFLNIYCDYE PGSEYNLDSI AQSCNLLELQ FTFPQLCHAE VQVGDQLETV FLLSGNDPAI HLYKENEGH QFEEQPVENL FPELTNLTSV VLWLDVHNFP GTSRRLSALG CQSGYVRVAH VDQRSREVLQ MWSVLQDGPI SRVIVFSLSA AKETKDRPLQ DEYSVLVASM LEPVVYRDL LNRGLEDQLL LPGSDQFDSV LCSLVTDVDL DGRPEVLVAT YGQELLCYKY RGPESGLPEA QHGFHLLWQR SFSSPLLAMA HVDLTGDGLQ ELAVVSLKGV HILQHSLIQA SELVLTRLRH

image



Background

Alternative Names KICSTOR complex protein kaptin (Actin-associated protein 2E4)

Product Name: Recombinant Human KPTN Protein
Catalog #: PEH9684



Note

For Research Use Only , Not for Diagnostic Use