


Product Name: Recombinant Human KIFAP3 Protein
Catalog #: PEH9685



Summary

Name	KIFA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92845
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HPSEKALIVH YEVEATILGE MGDPMGLGERK ECQKIIRLKS LNANTDITSL ARKVVEECKL IHPSKLNEVE QLLYYLQNRD DLSGKEKKE KSSKPKDPPP FEGMEIDEVA NINDMDEYIE LLYEDIPDKV RGSALILQLA RNPDNLEELL LNETALGALA

image



Background

Alternative Names Kinesin-associated protein 3 (KAP-3) (KAP3) (Smg GDS-associated protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KIFAP3 Protein
Catalog #: PEH9685

