

Product Name: Recombinant Human KIF17 Protein
Catalog #: PEH9687



Summary

Name	KIF17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2E2
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGQEAAGHLL GEQNYLPQEE PQEVPLQGLL GLQDPFAEVE AKLARLSSTV ARTDAPQADV PKVPVQVPAP TDLLEPSDAR PEAEAADDFP PRPEVDLASE VALEVRTAE PGVWLEAQAP VALVAQPEPL PATAGVKRES VGMEVAVLTD DPLPVVDQQQ VLARLQLEQ QVVGGEQAKN KDLKEKHRR KRYADERRKQ LVAALQNSDE DSGDWVLLNV YDSIQEEVRA KSKLLEKMQR KLRAAEVEIK DLQSEFQLEK IDYLATIRRQ ERDSMLLQQL

image



Background

Alternative Names Kinesin-like protein KIF17 (KIF3-related motor protein)

Product Name: Recombinant Human KIF17 Protein
Catalog #: PEH9687



Note

For Research Use Only , Not for Diagnostic Use