



Product Name: Recombinant Human KIF18A Protein
Catalog #: PEH9688



Summary

Name	KI18A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NI77
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LMPSPKGGQH TLKSPPSQSV QLNDLSKEL QPIVYTPEDC RKAQNPSTV TLMKPSSFTT SFQAISNIN SDNCLKMLCE VAIPHNRRKE CGQEDLDSTF TICEDIKSSK CKLPEQESLP NDNKDILQRL DPSSFSTKHS MPVPSMVPSY MAMTTAAKRK RKLTSSTSNS SLTADVNSGF AKRVRQDNSS EKHLQENKPT MEHKRNICKI

image



Background

Alternative Names Kinesin-like protein KIF18A (Marrow stromal KIF18A) (MS-KIF18A)

Note

Product Name: Recombinant Human KIF18A Protein
Catalog #: PEH9688



For Research Use Only , Not for Diagnostic Use