


Product Name: Recombinant Human LARP1B Protein
Catalog #: PEH9692



Summary

Name	LAR1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q659C4
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGTGVQVYPV EEALLKEYIK RQIEYYFSVE NLERDFFLRG KMDEQGFLPI SLIAGFQRVQ ALTTNLNLIL EALKDSTEVE IVDEKMRKKI EPEKWPIPGP PPRSVPPTDF SQLIDCPEFV PGQAFCSHTE SAPNSPRIGS PLSPKKNSET SILQAMSRGL STSLPDL DSE PWIEVKKRHQ PAPVKLRRESV SVPEGSLNQL CSSEEP EQEE LDFLFDEEIE QIGRKNTFTD WSDNDS DYEI DDQDLNKILI VTQTTPYVKK HPGGDRTGTH MSRAKITSEL AKVINDGLYY YEQDLWMEED ENKHTAIKQE VENFKKLNLI SKEQFENLTP ELPFEPNQEV PVAPSQSRQG GVQGV LHIPK KDLTDELAQK LFDVSEITSA AMVHSLPTAV

image



Background

Alternative Names	La-related protein 1B (La ribonucleoprotein domain family member 1B) (La
--------------------------	--

Product Name: Recombinant Human LARP1B Protein
Catalog #: PEH9692



ribonucleoprotein domain family member 2) (La-related protein 2)

Note

For Research Use Only , Not for Diagnostic Use