


Product Name: Recombinant Human LCMT1 Protein
Catalog #: PEH9694



Summary

Name	LCMT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9UIC8
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDLLPSKYFE VDFPMIVTRK LHSICKPPL SSPILELHSE DTLQMDGHIL DSKRYAVIGA DLRDLSELEE KLKCNMNTQ LPTLLIAECV LVYMTPEQSA NLLKWAANSF ERAMFINYEQ VNMGDRFGQI MIENLRRRQC DLAGVETCKS LESQKERLLS NGWETASAVD MMELYNRLPR AEVSRIESLE

image



Background

Alternative Names	Leucine carboxyl methyltransferase 1 (EC 2.1.1.233) (Protein-leucine O-methyltransferase) ([Phosphatase 2A protein]-leucine-carboxy methyltransferase 1)
--------------------------	--

Note

Product Name: Recombinant Human LCMT1 Protein
Catalog #: PEH9694



For Research Use Only , Not for Diagnostic Use