


Product Name: Recombinant Human L1RE1 Protein
Catalog #: PEH9699



Summary

Name	LORF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UN81
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELREEGFRRS NYSELREDIQ TKGKEVENFE KNLEECITRI TNTEKCLKEL MELKTKAREL REECRLRSR CDQLEERVSA MEDEMNMEMKR EGKFREKRIK RNEQSLQEIW

image



Background

Alternative Names	LINE-1 retrotransposable element ORF1 protein (L1ORF1p) (LINE retrotransposable element 1) (LINE1 retrotransposable element 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use