


Product Name: Recombinant Human MBLAC1 Protein
Catalog #: PEH9704



Summary

Name	MBLC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A4D2B0
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLLVPDYPYS VVLLQGYAE PEGVGDAVRA DGSVTLVLPQ TRGPASSHRE SPRSGGAEA ALEEAARGPI LVDTGGPWAR EALLGALAGQ GVAPGDVTLV VGTHGHSCHI GNLGLFPGAA LLVSHDFCLP GGRLPHGLG EGQPLRLGPG LEVWATPGHG QQRDVSVVVA GTALGTVVVA GDVFERDGE

image



Background

Alternative Names Metallo-beta-lactamase domain-containing protein 1 (EC 3.-.-.-)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MBLAC1 Protein
Catalog #: PEH9704

