


Product Name: Recombinant Human METTL27 Protein
Catalog #: PEH9706



Summary

Name	MET27_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N6F8
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVRARVRAAH GIPDLAQKLH FYDRWAPDYD QDVATLLYRA PRLAVDCLTQ ALPGPPHSAL ILDVACGTGL VAAELRAPGF LQLHGVDGSP GMLEQAQAPG LYQRLSLCTL GQEPLPSPEG TFDVAVLIVGA LSDGQVPCNA IPELVHTKPG GLVCLTTRTN SSNLQYKEAL EATLDRLEQA GMWEGLVAWP VDLRWTAGSW

image



Background

Alternative Names	Methyltransferase-like protein 27 (Williams-Beuren syndrome chromosomal region 27 protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human METTL27 Protein
Catalog #: PEH9706

