


Product Name: Recombinant Human MAP3K3 Protein
Catalog #: PEH9710



Summary

Name	M3K3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99759
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPHSGVSRQV RIKASQSAGD INTIYQPPEP RSRHLSVSSQ NPGRSSPPPG YVPERQQHIA RQGSYTSINS EGEFIPETSE QCMLDPLSSA ENSLSGSCQS LDRSADSPSF RKSRRMSRAQS FPDNRQEYSD RETQLYDKGV

image



Background

Alternative Names	Mitogen-activated protein kinase kinase kinase 3 (EC 2.7.11.25) (MAPK/ERK kinase kinase 3) (MEK kinase 3) (MEKK 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use