


Product Name: Recombinant Human MAD2L2 Protein
Catalog #: PEH9711



Summary

Name	MD2L2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UI95
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FGQVVADVLC EFLEVAVHLI LYVREVPVG IFQKRKKYNV PVQMSCHPEL NQYIQDTLHC VKPLEKNDV EKVVVVILDK EHRPVEKFVF EITQPPLLSI SSDSLSSHVE QLLRAFILKI SVCDAVLDHN PPGCTFTVLV HTREAATRNM EKIQVIKDFP WILADEQDVH MHDPRLIPLK TMTSDILKMQ LYVEERAHKG

image



Background

Alternative Names Mitotic spindle assembly checkpoint protein MAD2B (Mitotic arrest deficient 2-like protein 2) (MAD2-like protein 2) (REV7 homolog) (hREV7)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MAD2L2 Protein
Catalog #: PEH9711

