


Product Name: Recombinant Human ABHD10 Protein
Catalog #: PEH9717



Summary

Name	ABHDA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NUJ1
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IPGYLSYMNG TKALAIIEFC KSLGHACIRF DYSGVGSSDG NSEESTLGKW RKDVLSIIDD LADGPQILVG SSLGGWLMLH AAARPEKVV ALIGVATAAD TLVTKFNQLP VELKKEVEMK GVWSMPSKYS EEGVYNVQYS FIKEAEHHCL LHSPIPVNCP IRLHGMKDD

image



Background

Alternative Names	Palmitoyl-protein thioesterase ABHD10, mitochondrial (EC 3.1.2.22) (Acyl-protein thioesterase ABHD10) (Alpha/beta hydrolase domain-containing protein 10) (Abhydrolase domain-containing protein 10) (Mycophenolic acid acyl-glucuronide esterase, mitochondrial) (EC 3.1.1.93)
--------------------------	---

Product Name: Recombinant Human ABHD10 Protein
Catalog #: PEH9717



Note

For Research Use Only , Not for Diagnostic Use