


Product Name: Recombinant Human NAA15 Protein
Catalog #: PEH9721



Summary

Name	NAA15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXJ9
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IERHFIEITD DQDFH TYCM RKITLRSYVD LLKLEDVLRQ HPFYFKAARI AIEIYKLKLD NPLTDENKEH EADTANMSDK ELKCLRNKQR RAQKKAQIEE EKKNAEKEKQ QRNQKKKKDD

image



Background

Alternative Names	N-alpha-acetyltransferase 15, NatA auxiliary subunit (Gastric cancer antigen Ga19) (N-terminal acetyltransferase) (NMDA receptor-regulated protein 1) (Protein tubedown-1) (Tbdn100)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NAA15 Protein
Catalog #: PEH9721

