


**Product Name: Recombinant Human NECAB3 Protein**  
**Catalog #: PEH9726**



## Summary

<b>Name</b>	NECA3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96P71
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	19 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	FQDVFRRADK NDDGKLSFEE FQNYFADGVL SLGELQELFS GIDGHLTDNL ETEKLCDYFS EHLGVYRPVL AALES LNRAV LAAMDATKLE YERASKVDQF VTRFLLRETV SQLQALQSSL EGASDTLEAQ AHGWRSDAES VEAQSRLCGS

## image



## Background

**Alternative Names** N-terminal EF-hand calcium-binding protein 3 (Amyloid-beta A4 protein-binding family A member 2-binding protein) (Nek2-interacting protein 1) (Neuronal calcium-binding protein 3) (X11L-binding protein 51)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human NECAB3 Protein**  
**Catalog #: PEH9726**

---

