


Product Name: Recombinant Human METTL4 Protein
Catalog #: PEH9727



Summary

Name	METL4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N3J2
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NQGELDAMEY HTKIRELILD GSLLIQEGE KSGFLYPLFE KQDKGSKPIT LPLDACSLSE LCEMAKHLPS LNEMEHQTLQ LVEEDTSVTE QDLFLRVVEN NSSFTKVITL MGQKYLLPPK SSFLLSDISC MQPLLNYRKT FDVIVIDPPW

image



Background

Alternative Names	N(6)-adenine-specific methyltransferase METTL4 (Methyltransferase-like protein 4) (N(6)-adenine-specific DNA methyltransferase METTL4) (EC 2.1.1.72) (snRNA (2'-O-methyladenosine-N(6)-)-methyltransferase METTL4) (EC 2.1.1.-)
--------------------------	---

Note

Product Name: Recombinant Human METTL4 Protein
Catalog #: PEH9727



For Research Use Only , Not for Diagnostic Use