

Product Name: Recombinant Human NDUFV3 Protein
Catalog #: PEH9729



Summary

Name	NDUV3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P56181
Host	E.coli
Species	Human
Molecular Mass	9 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAESGKSEKG QPQNSKKQSP PKKPAPVPAE PFDNTTYKNL QHHDYSTYTF LDLNLELSKF RMPQPSSGRE SPRH

image



Background

Alternative Names	NADH dehydrogenase [ubiquinone] flavoprotein 3, mitochondrial (Complex I-9kD) (CI-9kD) (NADH-ubiquinone oxidoreductase 9 kDa subunit) (Renal carcinoma antigen NY-REN-4)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use