


Product Name: Recombinant Human NOXO1 Protein
Catalog #: PEH9730



Summary

Name	NOXO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFA2
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSRAAGRLSI HSLEAQLSRC LQPFCTQDTR DRPFQAQAE SLDVLLRHPS GWWLVENEDR QTAWFPAPYL EEAAPGQGRE GGPSLGSSGP QFCASRAYES SRADELSVPA GARVRVLETS DRGWWLCRYG DRAGLLPAVL LRPEGLGALL SGTGFRGGDD PAGEARGFPE PSQATAPPPT

image



Background

Alternative Names NADPH oxidase organizer 1 (NADPH oxidase regulatory protein) (Nox organizer 1) (Nox-organizing protein 1) (SH3 and PX domain-containing protein 5)

Note

Product Name: Recombinant Human NOX01 Protein
Catalog #: PEH9730



For Research Use Only , Not for Diagnostic Use