


Product Name: Recombinant Human WWP1 Protein
Catalog #: PEH9732



Summary

Name	WWP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0M0
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLERVKEQLK LLENKNGIA QTGELTVVLD GLVIEQENIT NCSSPTIEI QENGDALEN GEPSARTTAR LAEGTNGID NHVPTSTLVQ NSCCSYVVNG DNTPESSPSQV AARPKNTPAP KPLASEPADD TVNGESSSFA PTDNASVTGT

image



Background

Alternative Names	NEDD4-like E3 ubiquitin-protein ligase WWP1 (EC 2.3.2.26) (Atrophin-1-interacting protein 5) (AIP5) (HECT-type E3 ubiquitin transferase WWP1) (TGIF-interacting ubiquitin ligase 1) (Tiul1) (WW domain-containing protein 1)
--------------------------	--

Note

Product Name: Recombinant Human WWP1 Protein
Catalog #: PEH9732



For Research Use Only , Not for Diagnostic Use