


Product Name: Recombinant Human NUP214 Protein
Catalog #: PEH9737



Summary

Name	NU214_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35658
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVMAGIGEEI AHFQKELEEL KARTSKACFQ VGTSEEMKML RTESDDLHTF LLEIKETTES LHGDISSLKT TLLEGFAGVE EAREQNERNR DSGYLHLLYK RPLDPKSEAQ LQEIRRLHQY VKFAVQDVND

image



Background

Alternative Names	Nuclear pore complex protein Nup214 (214 kDa nucleoporin) (Nucleoporin Nup214) (Protein CAN)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use