



Summary

Name	NRIP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQ35
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRKETDMREA ASLRQRRMK QAVQFIHKDS ADLLPLDGLK KLGSSKDMQP HNILQRRLME TNLSKLRSGP RVPWASKTNK LNQAKSEGLK KSEEDDMILV SCQCAGKDVK ALVDTGCLYN LISLACVDRL GLKEHVKSHK HEGEKLSLPR HLKVVGQIEH LVITLGSRL DCPAAVVDDN EKNLSLGLQT LRSLKCIINL DKHRLIMGKT DKEEIPFVET VSLNEDNTSE

image



Background

Alternative Names Nuclear receptor-interacting protein 3 (Sarcoma antigen NY-SAR-105)

Note

Product Name: Recombinant Human NRIP3 Protein
Catalog #: PEH9739



For Research Use Only , Not for Diagnostic Use