


Product Name: Recombinant Human EXOG Protein
Catalog #: PEH9740



Summary

Name	EXOG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2C4
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQPDGSAEKA VLEQFGFPLT GTEARCYTNH ALSYDQAKRV PRWVLEHISK SKIMGDADRK HCKFKPDPNI PPTFSAFNED YVGSWWSRGGH MAPAGNNKFS SKAMAETFYL SNIVPQDFDN NSGYWNRIEM YCRELTERFE DVWVVSGPLT LPQTRGDGKK IVSYQVIGED NVAVPSHLYK VILARRSSVS TEPLALGAFV VPNEAIGFQP QLTEFQVSLQ DLEKLSGLVF FPHLDRTSDI RNICSVDTCK LLDFQEFTLY LSTRKIEGAR SVLRLEKIME NLKNAEIEPD DYFMSRYEKK LEELKAKEQS

image



Background

Alternative Names	Nuclease EXOG, mitochondrial (EC 3.1.30.-) (Endonuclease G-like 1) (Endo G-like 1)
--------------------------	--

Product Name: Recombinant Human EXOG Protein
Catalog #: PEH9740



Note

For Research Use Only , Not for Diagnostic Use