


Product Name: Recombinant Human DDX21 Protein
Catalog #: PEH9741



Summary

Name	DDX21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q9NR30
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TKKVTKNEEP SEEEIDAPKP KKMKEKEMN GETREKSPKL KNGFPHPEPD CNPSEAASEE SNSEIEQEIP VEQKEGAFSN FPISSEITKL LKGRGVTFLL PIQAKTFHHV YSGKDLIAQA RTGTGKTFSF AIPLIEKLHG ELQDRKRGRA PQVLVLAPTR ELANQVSKDF SDITKKLSVA CFYGGTPYGG QFERMRNGID ILVGTPGRIK DHIQNGKLDL TKLKHVVLDL VDQMLDMGFA DQVEEILSVA YKKDSEDNPQ TLLFSATCPH WVFNVAKKYM KSTYEQVDLI GKKTQKTAIT VEHLAIKCHW TQRAAVIGDV IRVYSGHQGR TIIFCETKKE AQELSQNDAI KQDAQSLHGD IPQKQREITL KGFRNGSFGV LVATNVAARG LDIPEVDLVI QSSPPKDVES YIHRSGRTGR AGRTGVCICF YQHKEEYQLV QVEQKAGIKF KRIGVPSATE

image



Background

Product Name: Recombinant Human DDX21 Protein
Catalog #: PEH9741



Alternative Names

Nucleolar RNA helicase 2 (EC 3.6.4.13) (DEAD box protein 21) (Gu-alpha)
(Nucleolar RNA helicase Gu) (Nucleolar RNA helicase II) (RH II/Gu)

Note

For Research Use Only , Not for Diagnostic Use