

Product Name: Recombinant Human OTUD6A Protein
Catalog #: PEH9745



Summary

Name	OTU6A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L8S5
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELQAQIRSLK NSVPKTDKTK RKQLLQDVAR MEAEMAQKHR QELEKFQDDS SIESVVEDLA KMNLENRPPR SSKAHRKRER MESEERERQE SIFQAEMSEH LAGFKREEEE KLAAILGARG

image



Background

Alternative Names OTU domain-containing protein 6A (EC 3.4.19.12) (DUBA-2)

Note

For Research Use Only , Not for Diagnostic Use