


Product Name: Recombinant Human PDZRN4 Protein
Catalog #: PEH9750



Summary

Name	PZRN4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-mSUMO & C-His
Accession #	Q6ZMN7
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EFLEELNLEM LEEEHNEAMQ PTANEVEQPK KQEEEEGTTD TATSSSNHE KDSGVGRTDE SLRNDESSEQ ENAAEDPNST SLKSKRDLGQ SQDTLGSVEL QYNESLVSGE YIDSDCIGNP DEDCERFRQL LELKCKIRNH GEYDLYSSS TIECNQGEQE GVEHELQLLN EELRNIELEC QNIMQAHLRQ KVTDQYGDIW TLHDGGFRNY NTSIDMQRGK LDDIMEHPEK SDKDSSSAYN

image



Background

Alternative Names PDZ domain-containing RING finger protein 4 (Ligand of Numb protein X 4) (SEMACAP3-like protein)

Note

Product Name: Recombinant Human PDZRN4 Protein
Catalog #: PEH9750



For Research Use Only , Not for Diagnostic Use